

CROSS CHAIR

Stretcher/Chair



Cross Chair is a device that adopts innovative solutions for the transport of patients to be used in rescue operations on road and ambulances.

The device can be used in chair or stretcher configuration, adapting to wider rescue scenarios. To reach more versatility and safety, there is a dedicated trolley.

Specific Characteristics

Configuration change takes less than 2 seconds
4-segment PVC Mattress with removable cover
Polyethylene board easy to clean
Backrest adjustable in 5 positions
Telescopic handles
Auto positioning rear handles
Suitable for use with Cross Chair Trolley
Light and robust
The side bars <i>Safe bar</i> contain the patient also in Chair configuration
Wheels \varnothing 150 mm to increase sliding action

Standard equipment

ST70002A	Black belt 2 pcs, metal buckle, reflex
ST43014	Cross Chair mattress

Accessories

IF01043B	IV holder TRACK 4
CR00081	Cross Chair trolley
ST00497B	DNA strap belt w/integrated re-winding system
ST00498B	DNA strap thoracic belt w/re-winding system

Technical data ¹

Height ^(CM)	1430 mm
Width ^(CM/SM)	600 mm
Depth with all handles closed ^(CM)	1000 mm
Depth with all handles opened ^(CM)	1365 mm
Length with handles closed ^(SM)	1970 mm
Length with handles opened ^(SM)	2235 mm
Side bar dimensions	600x200 mm
Weight	26 kg
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

CROSS CHAIR
STRETCHER/CHAIR YELLOW/B
CR00080C



CND Classification

V0899

¹CM:Chair mode – SM: Stretcher mode

Rev.0 (28/10/2014)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it